

Canna

Cannulated Screw System



Features:

1. Complete ranges of cannulated screw can mathvarious clinical needs.



2. Stardrive head design prevents screw slippage.

3. Self-Tapping screw shorten surgery time.



4. Different thread length accommodates different fragment compression needs.

5. Variable thread pitch achieve fragment compression upon screw insertion.



6. Full Thread Compression Screws create compression and secure stability

7. Our products are manufactured using one or more of the following materials:commercially pure titanium(IOS 5832-2,ASTM F67),Ti6Al4V alloy(ISO 5832-3,ASTM F1472).

Standard Cannulated Screw-Φ2.4

Product Code	Specification		
	length	thread	
52702410	10		
52702411	11		
52702412	12		
52702413	13		
52702414	14		
52702415	15		
52702416	16		
52702417	17		
52702418	18		
52702419	19		
52702420	20		short
52702421	21		
52702422	22		
52702423	23		
52702424	24		
52702425	25		
52702426	26		
52702427	27		
52702428	28		
52702429	29		
52702430	30		
52712417	17		
52712418	18		
52712419	19		
52712420	20		
52712421	21		
52712422	22		
52712423	23		
52712424	24		long
52712425	25		
52712426	26		
52712427	27		
52712428	28		
52712429	29		
52712430	30		



Standard Cannulated Screw-Φ3.0

Product Code	Specification		
	length	thread	
52703008	8		
52703009	9		
52703010	10		
52703011	11		
52703012	12		
52703013	13		
52703014	14		
52703015	15		
52703016	16		
52703017	17		
52703018	18		
52703019	19		
52703020	20		
52703021	21		short
52703022	22		
52703023	23		
52703024	24		
52703025	25		
52703026	26		
52703027	27		
52703028	28		
52703029	29		
52703030	30		
52703032	32		
52703034	34		
52703036	36		
52703038	38		
52703040	40		
52713014	14		
52713015	15		
52713016	16		
52713017	17		
52713018	18		
52713019	19		
52713020	20		
52713021	21		
52713022	22		
52713023	23		
52713024	24		long
52713025	25		
52713026	26		
52713027	27		
52713028	28		
52713029	29		
52713030	30		
52713032	32		
52713034	34		
52713036	36		
52713038	38		
52713040	40		



Standard Cannulated Screw-Φ4.0

Product Code	Specification		
	length	thread	
52704010	10		
52704012	12		
52704014	14		
52704016	16		
52704018	18		
52704020	20		
52704022	22		
52704024	24		
52704026	26		
52704028	28		
52704030	30		
52704032	32		
52704034	34		
52704036	36		
52704038	38		short
52704040	40		
52704042	42		
52704044	44		
52704046	46		
52704048	48		
52704050	50		
52704052	52		
52704054	54		
52704056	56		
52704058	58		
52704060	60		
52704064	64		
52704068	68		
52704072	72		
52714016	16		
52714018	18		
52714020	20		
52714022	22		
52714024	24		
52714026	26		
52714028	28		
52714030	30		
52714032	32		
52714034	34		
52714036	36		
52714038	38		
52714040	40		long
52714042	42		
52714044	44		
52714046	46		
52714048	48		
52714050	50		
52714052	52		
52714054	54		
52714056	56		
52714058	58		
52714060	60		
52714064	64		
52714068	68		
52714072	72		




Headless Compression Screw-Φ4.5

Product Code	Specification	
	length	thread
52804526	26	
52804528	28	
52804530	30	
52804532	32	
52804534	34	
52804536	36	
52804538	38	
52804540	40	
52804542	42	
52804544	44	
52804546	46	
52804548	48	
52804550	50	
52804552	52	
52804554	54	
52804556	56	
52804558	58	
52804560	60	
52804564	64	
52804568	68	
52804572	72	
52804576	76	
52804580	80	
52814530	30	
52814532	32	
52814534	34	
52814536	36	
52814538	38	
52814540	40	
52814542	42	
52814544	44	
52814546	46	
52814548	48	
52814550	50	
52814552	52	
52814554	54	
52814556	56	
52814558	58	
52814560	60	
52814564	64	
52814568	68	
52814572	72	
52814576	76	
52814580	80	

short

long



Headless Compression Screw-Φ5.5

Product Code	Specification	
	length	thread
52805525	25	
52805530	30	
52805535	35	
52806540	40	
52806545	45	
52806550	50	
52806555	55	
52806560	60	
52806565	65	
52806570	70	
52815540	40	
52815545	45	
52815550	50	
52815555	55	
52815560	60	
52815565	65	
52815570	70	
52815575	75	
52815580	80	
52815585	85	
52815590	90	
52815595	95	
52815500	100	

short

long



Headless Compression Screw-Φ6.5

Product Code	Specification	
	length	thread
52806530	30	
52806535	35	
52806540	40	
52806545	45	
52806550	50	
52806555	55	
52806560	60	
52806565	65	
52806570	70	
52806575	75	
52806580	80	
52806585	85	
52806590	90	
52806595	95	
52806500	100	
52806505	105	
52806510	110	
52806515	115	
52806520	120	
52816550	50	
52816555	55	
52816560	60	
52816565	65	
52816570	70	
52816575	75	
52816580	80	
52816585	85	
52816590	90	
52816595	95	
52816500	100	
52816505	105	
52816510	110	
52816515	115	
52816520	120	

short

long



Full Thread Compression Screw-Φ2.5

Product Code	Specification
	length
52902508	8
52902510	10
52902512	12
52902514	14
52902516	16
52902518	18
52902520	20
52902522	22
52902524	24
52902526	26



Full Thread Compression Screw-Φ4.0

Product Code	Specification
	length
52904016	16
52904018	18
52904020	20
52904022	22
52904024	24
52904026	26
52904028	28
52904030	30
52904032	32
52904034	34
52904036	36
52904038	38
52904040	40
52904042	42
52904044	44
52904046	46
52904048	48
52904050	50



Standard Washers

Product Code	Specification
52740045	Φ4.5
52740065	Φ6.5
52740070	Φ7.0
52740100	Φ10
52740130	Φ13
52740160	Φ16



Full Thread Compression Screw-Φ3.5

Product Code	Specification
	length
52903516	16
52903518	18
52903520	20
52903522	22
52903524	24
52903526	26
52903528	28
52903530	30
52903532	32
52903534	34
52903536	36
52903538	38
52903540	40



Full Thread Compression Screw-Φ5.2

Product Code	Specification
	length
52905225	25
52905230	30
52905235	35
52905240	40
52905245	45
52905250	50
52905255	55
52905260	60
52905265	65
52905270	70
52905275	75
52905280	80



Deflection Washers

Product Code	Specification
52750130	Φ13

