

Ramic

- Intramedullary Nail Systems
- Femoral Nail



Features

1. Anatomical Design

--Optimal fit in bone



--Easy insertion and extraction

--Flat lateral side to reduce the damage on lateral cortex



2. Multiple inter-locking option

--Standard locking

--Recon locking

--Distal multi-optional inter-locking design to increase stability

--Distal dynamic locking option

3. Cannulated Stardrive End Cap makes easy insertion and extraction.

4. Innovative and friendly-used instrumentation.

5. Our products are manufactured using one or more of the following materials: commercially pure titanium (ISO 5832-2, ASTM F67), Ti6Al4V alloy (ISO 5832-3, ASTM F1472).

Indications

1. Fractures of femoral shaft
2. Combined femoral shaft and neck fractures (ipsilateral)
3. Sub-trochanteric fractures

Fomoral Nail

Product Code Left	Product Code Right	Diameter	Length
55009034	55019034	φ9.0	340
55009036	55019036	φ9.0	360
55009038	55019038	φ9.0	380
55009040	55019040	φ9.0	400
55009042	55019042	φ9.0	420
55001034	55011034	φ10	340
55001036	55011036	φ10	360
55001038	55011038	φ10	380
55001040	55011040	φ10	400
55001042	55011042	φ10	420
55001134	55011134	φ11	340
55001136	55011136	φ11	360
55001138	55011138	φ11	380
55001140	55011140	φ11	400
55001142	55011142	φ11	420
55001234	55011234	φ12	340
55001236	55011236	φ12	360
55001238	55011238	φ12	380
55001240	55011240	φ12	400
55001242	55011242	φ12	420

Locking screw-A

Product Code	Diameter	Length
55030026	φ5.0	26
55030028	φ5.0	28
55030030	φ5.0	30
55030032	φ5.0	32
55030034	φ5.0	34
55030036	φ5.0	36
55030038	φ5.0	38
55030040	φ5.0	40
55030042	φ5.0	42
55030044	φ5.0	44
55030046	φ5.0	46
55030048	φ5.0	48
55030050	φ5.0	50
55030052	φ5.0	52
55030054	φ5.0	54
55030056	φ5.0	56
55030058	φ5.0	58
55030060	φ5.0	60
55030062	φ5.0	62
55030064	φ5.0	64
55030066	φ5.0	66
55030068	φ5.0	68
55030070	φ5.0	70
55030072	φ5.0	72
55030074	φ5.0	74
55030076	φ5.0	76
55030078	φ5.0	78
55030080	φ5.0	80

Fomoral Hip Screws

Product Code	Diameter	Length
55020060	φ6.5	60
55020065	φ6.5	65
55020070	φ6.5	70
55020075	φ6.5	75
55020080	φ6.5	80
55020085	φ6.5	85
55020090	φ6.5	90
55020095	φ6.5	95
55020100	φ6.5	100
55020105	φ6.5	105
55020110	φ6.5	110
55020115	φ6.5	115
55020120	φ6.5	120
55020125	φ6.5	125
55020130	φ6.5	130

Fomoral End Caps

Product Code	Length
55070000	0
55070005	5
55070010	10
55070015	15
55070020	20

Locking screw-B

Product Code	Diameter	Length
55030026	φ5.0	26
55030028	φ5.0	28
55030030	φ5.0	30
55030032	φ5.0	32
55030034	φ5.0	34
55030036	φ5.0	36
55030038	φ5.0	38
55030040	φ5.0	40
55030042	φ5.0	42
55030044	φ5.0	44
55030046	φ5.0	46
55030048	φ5.0	48
55030050	φ5.0	50